

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Fax: 440/337-1180

web: www.mpbio.com

Certificate of Analysis

Product Description: Peroxidase-Conjugated Rabbit IgG Fraction To Bovine Albumin
Catalog Number: 55285
Lot: S3597

Formula: N/A
CAS #: N/A
Physical Description: Brown Solid

Formula Weight: N/A
Storage: 2 - 8°C

Test	Specification	Result
Identity Test	Passes	Passes
Purity	≥ 95%	≥95% (SDS-Page)
Concentration	Report Result	30 mg/ml
Formulation	PBS (no preservative)	PBS (no preservative)

Specificity:

- Bovine Albumin Purified Protein: one band observed
- Normal Bovine Serum: one band observed
- Bovine IgG and Bovine Transferrin: negative for both

ELISA: Expected range of 1:32K - 1:64K

Store lyophilized product at 2-8°C. For long-term storage, aliquot and freeze reconstituted product at -20°C or lower.

Recommended Retest Date (MM/DD/YYYY): 02/12/2023

02/14/2020 - John Huang, PhD
 MP Biomedicals, LLC.
 Quality Control Manager

This is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
 Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)