

**MP Biomedicals, LLC**29525 Fountain Parkway  
Solon, Ohio 44139Telephone: 440/337-1200  
Toll Free: 800/854-0530Fax: 440/337-1180  
web: www.mpbio.com

## Certificate of Analysis

**Product Description:** Anti-human Synaptophysin  
**Catalog Number:** 10706  
**Lot:** S4067

**Formula:** N/A  
**CAS #:** N/A  
**Physical Description:** Clear, colorless liquid

**Formula Weight:** N/A  
**Storage:** 2 - 8°C

Test	Specification	Result
Host	Rabbit	Rabbit
Form	Liquid	Liquid
Antigenic Specificity	Synaptophysin	Synaptophysin
Immunogen	Synthetic human synaptophysin	Synthetic human synaptophysin

Concentration: 3.0 mg/ml

Buffer: 10 mM PO<sub>4</sub>, 150mM NaCl, pH7.4 containing 1% BSA as stabilizer and 0.05% Sodium Azide as a preservative

Recommended Retest Date (MM/DD/YYYY): 04/20/2021

A handwritten signature in black ink that reads 'John Huang'.

04/23/2020 - John Huang, PhD  
 MP Biomedicals, LLC.  
 Quality Control Manager

This is an electronically generated document  
<mailto:biotech@mpbio.com>  
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site  
 Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)