

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: Mouse, Anti-beta-Galactosidase, Monoclonal Antibody, Purified
Catalog Number: 63365
Lot: S4249

Formula: N/A
CAS #: N/A
Physical Description: Frozen Liquid

Formula Weight: N/A
Storage: -20°C

Test	Specification	Result
Identity Test	Passes	Passes
Source	Mouse	Mouse
Form	Liquid	Liquid
Isotype	IgG2a(k)	IgG2a(k)
Antigenic Specificity	beta- Galactosidase	beta- Galactosidase
Clone	6A-11	6A-11
Concentration	For Information Only	2.0-2.5 mg/ml in a solution containing 100mM Tris-HCl (pH 8.0), 150 mM NaCl, 0.02% NaN ₃ .

Purity: ≥95% (SDS-Polyacrylamide gel electrophoresis)
Recommended Retest Date (MM/DD/YYYY): 06/02/2023

06/04/2020 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
mailto:biotech@mpbio.com
http://www.mpbio.com

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)