

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530Fax: 440/337-1180
web: www.mpbio.com

Certificate of Analysis

Product Description: Peroxidase-Conjugated Rabbit IgG Fraction To Bovine Albumin Catalog Number: 55285 Lot: S4587
--

Formula: N/A CAS #: N/A Physical Description: Brown Solid
--

Formula Weight: N/A Storage: 2 - 8°C

Test	Specification	Result
Identity Test	Passes	Passes
Purity	≥ 95%	≥ 95% (SDS-Page)
Concentration	Report Result	30 mg/ml
Formulation	PBS (no preservative)	PBS (no preservative)

Specificity:

- Bovine Albumin Purified Protein: one band observed
- Normal Bovine Serum: one band observed
- Bovine IgG and Bovine Transferrin: negative for both

ELISA: Expected range of 1:64K - 1:128K

Store lyophilized product at 2-8°C. For long-term storage, aliquot and freeze reconstituted product at -20°C or lower.

Recommended Retest Date (MM/DD/YYYY): 08/06/2023

08/10/2020 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control ManagerThis is an electronically generated document
mailto:biotech@mpbio.com
http://www.mpbio.comOnline Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)