

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: Human anti-SARS-CoV-2 IgG ELISA Kit
Catalog Number: 440100
Lot: S5390**Formula:** N/A
CAS #: N/A
Physical Description: Kit**Formula Weight:** N/A
Storage: 2 - 8°C

Component	Cat#	UN#, Packing Group, Classification	Volume
Human anti SARS-CoV2(N) IgG Coated Micro Plate	E001		8 x 12
Human anti SARS-CoV2(N) IgG Standard	E002		500 ul
Sample Dilution Buffer	E039		20 mL
25 x Wash Buffer	E038		30 mL
TMB Substrate	E024		5 mL
Stop Solution	E026		5 mL
Antibody Dilution Buffer	E040		5 mL
HRP-Anti Human IgG Antibody	E003		60 ul

Test	Specification	Result
Identity Test	Passes	Passes

Recommended Retest Date (MM/DD/YYYY): 06/14/2021

A handwritten signature in black ink, appearing to read 'J Huang'.

01/12/2021 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control ManagerThis is an electronically generated document
mailto:biotech@mpbio.com
http://www.mpbio.comOnline Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)