

**MP Biomedicals, LLC**29525 Fountain Parkway  
Solon, Ohio 44139Telephone: 440/337-1200  
Toll Free: 800/854-0530

## Certificate of Analysis

**Product Description:** Antigen affinity gel rabbit IgG whole molecule  
**Catalog Number:** 55961  
**Lot:** S5724

**Formula:** N/A  
**CAS #:** N/A  
**Physical Description:** Cloudy White suspension/ solution

**Formula Weight:** N/A  
**Storage:** 2 - 8°C

Test	Specification	Result
Identity Test	Passes	Passes

Buffer: 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.5

Preservative: 0.05% sodium azide

Antibody Concentration: 1.0 mg/mL

Purity: Based on immunoelectrophoresis at an antigen concentration of 20 mg/mL, the pattern of precipitation against anti-rabbit whole serum is the same as that against anti-rabbit IgG.

Do not freeze.

For research or further manufacturing use only. Not intended for human or animal therapeutic or diagnostic procedures.

Recommended Retest Date (MM/DD/YYYY): 06/20/2022

03/03/2021 - John Huang, PhD  
MP Biomedicals, LLC.  
Quality Control Manager

This is an electronically generated document  
<mailto:biotech@mpbio.com>  
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site  
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)