

## Certificate of Analysis

**Product Description:** Replacement O-rings for Metal Lysing Tube, 2 mL  
**Catalog Number:** 6990100  
**Lot:** S5789

**Formula:** N/A  
**CAS #:** N/A  
**Physical Description:** Replacement O-rings for Metal Lysing Tube, 2 mL

**Formula Weight:** N/A  
**Storage:** 15 - 30°C

Test	Specification	Result
Identity Test	Passes	Passes

Recommended Retest Date (MM/DD/YYYY): 05/04/2026



05/05/2021 - John Huang, PhD  
MP Biomedicals, LLC.  
Quality Control Manager

This is an electronically generated document  
<mailto:biotech@mpbio.com>  
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site  
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)