

Certificate of Analysis

Product Description: Rabbit, Anti-Toxoplasma Gondii C56
Catalog Number: 68813
Lot: S6661

Formula: N/A
CAS #: N/A
Physical Description: Clear Colorless Liquid

Formula Weight: N/A
Storage: 2 - 8°C, do not freeze

Test	Specification	Result
Identity Test	Passes	Passes
Host	Rabbit	Rabbit
Form	Liquid	Liquid
Antigenic Specificity	Toxoplasma gondii	Toxoplasma gondii
Immunogen	Toxoplasma gondii strain C56	Toxoplasma gondii strain C56

Intended Use: This antibody is currently available for research use only. This antibody is designed for the specific localization of *Toxoplasma gondii* in formalin-fixed, paraffin-embedded tissue sections.

Reagents Provided: Pooled filtered rabbit antisera diluted in phosphate buffered saline, pH 7.6, containing 1% BSA and 0.09% sodium azide.

Storage and Handling: Antibodies should be stored at 2 - 8°C without further dilution. Fresh dilutions, if required, should be made prior to use and are stable for up to one day at 15 - 30°C. Unused portions of antibody preparations should be discarded after one day to minimize microbial contamination and increase in nonspecific staining. Diluted antibody preparations can be refrigerated or frozen for extended shelf life.

Recommended Retest Date (MM/DD/YYYY): 07/31/2023



07/21/2021 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
mailto:biotech@mpbio.com
http://www.mpbio.com

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)