

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: TRIS-Glycine-Native Running Buffer (10X)
Catalog Number: 822521
Lot: S7208

Formula: N/A
CAS #: N/A
Physical Description: Clear Colorless Liquid

Formula Weight: N/A
Storage: 15°C - 30°C

Test	Specification	Result
Identity Test	Passes	Passes
Expiration Date	mm/dd/yyyy	09/30/2022

Tris-Base: 250mM
Glycine: 1.92M
pH(1X)=8.5 ± 0.15

10/06/2021 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)