



Certificate of Analysis

Product Description: 24 x 2 ml CoolPrep Sample Holder
Catalog Number: 6002528
Lot: S7493

Formula: N/A
CAS #: N/A
Physical Description: 24 x 2 ml Holder

Formula Weight: N/A
Storage: 15 - 30°C

Test	Specification	Result
Identity Test	Passes	Passes

Recommended Retest Date (MM/DD/YYYY): 11/04/2031

A handwritten signature in black ink, appearing to read 'J Huang'.

11/08/2021 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)