

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: SafTest analyzer replacement bulb (light source)
Catalog Number: WS1024
Lot: S8072**Formula:** N/A
CAS #: N/A
Physical Description: 2 pin cable attached to flat board with LED light**Formula Weight:** N/A
Storage: 15 - 30°C

Test	Specification	Result
Identity Test	Passes	Passes

Recommended Retest Date (MM/DD/YYYY): 03/07/2027

03/08/2022 - Fengqi An
MP Biomedicals, LLC.
Life Science QC SupervisorThis is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)