

## Certificate of Analysis

**Product Description:** Mouse anti-HA Tag mAb  
**Catalog Number:** L100032  
**Lot:** S9163

**Formula:** N/A  
**CAS #:** N/A  
**Physical Description:** Liquid

**Formula Weight:** N/A  
**Storage:** -20°C

Test	Specification	Result
Identity Test	Passes	Passes
Host	Mouse	Mouse
Form	Liquid	Liquid
Clone	Monoclonal Antibody	Monoclonal Antibody
Conjugate	Unconjugated	Unconjugated
Antigenic Specificity	HA Tag	HA Tag
Concentration	Report Result	1.04 mg/mL

**DESCRIPTION:** Immunogen is a synthetic peptide corresponding to HA tag, this affinity-purified antibody is prepared by solid phase absorption. The resulting antibody is highly pure and offers high specificity and potency against the tag antigen. The antibody is dissolved in PBS buffer with 0.02% sodium azide, 1% BSA, 50% glycerol, pH 7.3.

**USES:** This monoclonal antibody has been experimentally validated in the following applications by offering high specificity and potency against the tag antigen.

**WB:** 1:2000 - 1:5000; **IF:** 1:50 - 1:500; **IP:** 1:20 - 1:50.

**STORAGE:** This antibody can be stored at 2-8°C for one month without detectable loss of activity. This antibody is also stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

**WARNING AND CAUTION:** This product is intended FOR RESEARCH USE ONLY, not for use in diagnostic or therapeutic procedures involving humans.

**Recommended Retest Date (MM/DD/YYYY):** 07/27/2025

A handwritten signature in black ink on a light gray grid background. The signature appears to be 'F. An'.

08/10/2022 - *Fengqi An*  
MP Biomedicals, LLC.  
Life Science QC Supervisor

This is an electronically generated document  
<mailto:biotech@mpbio.com>  
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site  
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)