

**MP Biomedicals, LLC**29525 Fountain Parkway  
Solon, Ohio 44139Telephone: 440/337-1200  
Toll Free: 800/854-0530

## Certificate of Analysis

**Product Description:** Pectolyase  
**Catalog Number:** 151804  
**Lot:** SR01725**Formula:** N/A  
**CAS #:** N/A  
**Physical Description:** Tan Powder**Formula Weight:** N/A  
**Storage:** 2 - 8°C

Test	Specification	Result
Identity Test	Passes	Passes
Source	<i>Aspergillus japonicus</i>	<i>Aspergillus japonicus</i>
Activity	>1 u/mg solid	>1,000 u/g (pectin lyase)

One unit will liberate 1.0  $\mu$ mole of galacturonic acid from polygalacturonic acid per minute at pH 5.5 at 25°C.

Form: Lyophilized

Pectinase types: Edopolygalacturonase and Endo-pectin lyase

Recommended Retest Date (MM/DD/YYYY): 08/06/2024



10/01/2019 - John Huang, PhD  
MP Biomedicals, LLC.  
Quality Control Manager

This is an electronically generated document  
mailto:biotech@mpbio.com  
http://www.mpbio.com

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site  
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)