

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530Fax: 440/337-1180
web: www.mpbio.com

Certificate of Analysis

Product Description: Anti-human Synaptophysin
Catalog Number: 10706
Lot: SR03529

Formula: N/A
CAS #: N/A
Physical Description: Clear, colorless liquid

Formula Weight: N/A
Storage: 2 - 8°C

Test	Specification	Result
Host	Rabbit	Rabbit
Form	Liquid	Liquid
Antigenic Specificity	Synaptophysin	Synaptophysin
Immunogen	Synthetic human synaptophysin	Synthetic human synaptophysin

Concentration: 3.0 mg/ml

Buffer: 10 mM PO₄, 150mM NaCl, pH7.4 containing 1% BSA as stabilizer and 0.05% Sodium Azide as a preservative

Recommended Retest Date (MM/DD/YYYY): 06/20/2022

A handwritten signature in black ink that reads 'John Huang'.

12/21/2020 - John Huang, PhD
 MP Biomedicals, LLC.
 Quality Control Manager

This is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
 Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)