

MP Biomedicals, LLC29525 Fountain Parkway
Solon, Ohio 44139Telephone: 440/337-1200
Toll Free: 800/854-0530

Certificate of Analysis

Product Description: Aprotinin Solution, 200,000 KIU/ml
Catalog Number: 180563
Lot: SR04215**Formula:** N/A
CAS #: 9087-70-1
Physical Description: Clear, Very Light Brown Solution**Formula Weight:** N/A
Storage: 2 - 8°C

Test	Specification	Result
Identity Test	Passes	Passes
Source	Bovine Lung	Bovine Lung
Activity	>200,000 KIU/ml	214,200 KIU/ml
Sodium Chloride	~9 mg/ml	9 mg/ml
Benzyl Alcohol	~9 mg/ml	10 mg/ml
pH (For Information Only)	Report Result	5.68
Expiration Date	Month/Day/Year	08/03/2023

05/28/2021 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control ManagerThis is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)