

Certificate of Analysis

Product Description: Rabbit, Anti-Bovine Albumin
Catalog Number: 65111
Lot: SR04261

Formula: N/A
CAS #: N/A
Physical Description: Lyophilized

Formula Weight: N/A
Storage: 2 - 8°C

Test	Specification	Result
Identity Test	Passes	Passes
Host	Rabbit	Rabbit
Form	Lyophilized	Lyophilized
Antigenic Specificity	Bovine Albumin	Bovine Albumin
Concentration	For Information Only	39 mg/mL

Procedure: Rabbits were immunized with highly purified bovine serum albumin (BSA) and the resulting antiserum was collected. The antiserum was delipidated, solid phase adsorbed and 0.45um filtered. Total protein measurement was determined by refractometer. Specific antibody content is determined by immunoaffinity purification off an antigen containing immunosorbent. Specifications: This antiserum reacts with BSA as determined by IEP and Ouchterlony techniques. This antiserum may cross react with albumin from other species. Optimal working dilutions should be determined experimentally by the investigator.

For research and manufacturing only; not for use in diagnostic applications. For in vitro use only, not for bovine or animal consumption.

Recommended Retest Date (MM/DD/YYYY): 12/17/2022



06/17/2021 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
mailto:biotech@mpbio.com
http://www.mpbio.com

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)