



**MP Biomedicals, LLC**

29525 Fountain Parkway  
Solon, Ohio 44139

Telephone:  
440/337 · 1200  
Toll Free: 800/854 · 0530  
Fax: 440/337 · 1180

mailto:biotech@mpbio.com  
web: www.mpbio.com

# cappel

**Goat Anti-Chicken Serum (Forensic)**  
**2ML**  
**Catalog #: 56000**  
**Lot #: U1123173250**

Total Protein: 63.8 mg/ml

**DESCRIPTION:**

Product is the lyophilized powder of goat anti-chicken serum (forensic) and buffer salts.

The antiserum is prepared by immunizing goats with chicken serum; antiserum is delipidated and made specific by immunoabsorption with soluble antigens or solid-phase immunoabsorption where appropriate. The resulting antiserum is dialyzed into 0.02M sodium phosphate, 0.14M sodium chloride, pH 7.3, with 0.05% sodium azide, adjusted to standard titer, filtered through a 0.22 µm filter, vialled and lyophilized.

**USES:**

This product is intended for use as a precipitating antibody. It gives sharp precipitin lines in immunoelectrophoresis, capillary tube, and immunodiffusion techniques.

**SPECIFICATIONS:**

The total protein is measured using the Biuret procedure, with bovine albumin as the standard. The product shows the following reactions at the stated chicken serum dilutions and the following reaction to turkey serum because of species-homology.

Chicken Serum Dilution	Microtube Precipitin Assay	
	15 Minutes	30 Minutes
1:3	++++	++++
1:10	++	+++
1:100	+	++
1:1000	O	+

The product has been tested at 60 mg/ml and found negative against the following whole sera at 20-25 mg/ml:

- |        |            |        |        |
|--------|------------|--------|--------|
| Bear   | Dog        | Horse  | Rabbit |
| Bovine | Goat       | Human  | Rat    |
| Cat    | Guinea Pig | Monkey | Sheep  |
| Deer   | Hamster    | Mouse  | Swine  |

**EXPIRATION DATE:**

07/31/2028

**RECONSTITUTION:**

Reconstitute product with 2.0 ml of deionized or distilled water. Gentle swirling may be used to speed rehydration. Avoid vigorous shaking of the reconstituted material. If a small amount of precipitate is present, or if preparation appears cloudy, clarify by centrifugation or filtration.

**STORAGE:**

Store product at 2-8°C. Store reconstituted material no longer than 4 weeks at 2-8°C. For long-term storage, aliquot and freeze reconstituted product at -20°C or lower.

**WARNING AND CAUTION:**

This product is intended **FOR RESEARCH USE ONLY**, not for use in diagnostic or therapeutic procedures involving humans or animals.

This product contains sodium azide. To prevent formation of toxic vapors, do not mix with strongly acidic solutions. To prevent formation of potentially explosive metallic azides in metal plumbing, always wash into drain with copious quantities of water.

**RELATED PRODUCT ORDERING INFORMATION:**

<u>Catalog #</u>	<u>Description</u>	<u>Size</u>
55830	ID DISCS PATTERN-A WITH INDICATOR	20EA
55831	ID DISCS PATTERN-A WITHOUT INDICATOR	20EA

**For whole sera and other forensic antisera, please see our catalog.**



**Approved by:** *Fengqi An*  
Quality Control Manager

**TO PLACE AN ORDER IN THE U.S., PLEASE CALL 800-854-0530  
OR TO OBTAIN TECHNICAL ASSISTANCE IN THE U.S., PLEASE CALL 800-279-5490  
OUTSIDE THE U.S., CONTACT YOUR LOCAL DISTRIBUTOR.**